

MRI BRAIN STROKE PROTOCOL

Patient Position Supine

Localizer scout

Plane	Sequence	Fat Sat	Mode	SLT / SP	FOV	Notes
Axial	DWI			5/1.5	24 cm	
Sagittal	T1			5/1.5	24 cm	
Axial	T1			5/1.5	24 cm	
Axial	T2			5/1.5	24 cm	
Axial	T2 FLAIR			5/1.5	24 cm	
Axial	**GRE			5/1.5	24 cm	**3D SWAN preferred; GRE if time constraints/pt motion

**Perform your scanner equivalent

axial whole head



3D TOF slab

