Reviewed: Jan 2022 Updated: Oct 2020



MRA NECK WITHOUT & WITH CONTRAST

Indications	Intracranial arterial stenosis, occulsion, abnormality		
IV site	IV site in RIGHT arm		
Localizer	3 plane scout		
Notes	Precontrast MRA necessary before post contrast performed		

PLANE	SEQUENCE	FAT SAT	SLT/SP	FOV	NOTES
Axial	T1	Υ	3-5 mm	20 cm	thyroid gland to skull base
Axial	2D and/or 3D TOF		< 1.5 or VOI	24 cm	cover from thoracic inlet to skull base
Axial	3D TOF		VOI	24 cm	minimum coverage from aortic arch to floor of sella. Subclavian arteries 2 cm distal to the vertebral arteries.

^{**}Do vendor appropriate timing sequence prior to CE-MRA

Interval ≤0.75 on recons

Post processing

- · MIP right half of neck in R to L rotation
- MIP Left half of neck in R to L rotation
- MIP entire neck in R to L rotation
- MPR in sagittal and coronal plane. Send source images



