

MRI PELVIS WITHOUT - SPORTS HERNIA

Patient Position	Supine-arms out of FOV
FOV	Fits to patient/ body part (~350 mm)
Localizer	3 plane scout localize
Scan Range	Middle of lower kidney through symphysis pubis
Breathing	Scan on Expiration

PLANE	SEQUENCE	FAT SAT	MODE	SLT/SP	FOV	NOTES
Axial	T2 SPAIR		SPAIR	3 / 1	both hips	1 cm above acetabulum thru symphysis pubis
Axial OBL	PD (TE=20msec)			3 / 1	20 cm	parallel anterior iliac crest. From crest thru hip joint.
Axial OBL	T2 SPAIR		SPAIR	3/1	20 cm	parallel anterior iliac crest. From crest thru hip joint.
Sagittal	T2 SPAIR		SPAIR	3 / 1	20 cm	Ischial tuberosity to ischial tuberosity (NOT skin to skin)
Coronal	T1		None	3 / 1	both hips	1cm anterior to symphysis pubis thru both ischial tuberosities
Coronal	T2 STIR (TE= 50-60)		STIR	3 / 1	both hips	1cm anterior to symphysis pubis thru both ischial tuberosities

See next page for images

MRI PELVIS WITHOUT - SPORTS HERNIA

