

## MRI STERNOCLAVICULAR JOINTS without contrast (BILATERAL)

<b>Indications</b>	<b>Degenerative or inflammatory conditions, congenital anomalies, benign lesions.</b>
FOV	Include both joints and medial 1/3 of clavicles.
Position	PRONE - if possible with torso or large flex coil.
Localizer	3 plane scout localizer
Note	Always use marker at the site of most pain, injury, or swelling.

<b>PLANE</b>	<b>SEQUENCE</b>	<b>FAT SAT</b>	<b>FOV</b>	<b>NOTES</b>
Coronal	STIR			
Coronal	T1	N		
Axial	PD	Y		
Axial	T1	N		
Sagittal	T2	Y		